

Prijenosni stimulatori za ublažavanje boli

Za brzo i učinkovito ublažavanje
boli ili otklanjanja boli

- Elektroterapija
- Terapija ultrazvukom
- Kombinacija elektroterapije i ultrazvučne terapije



Medio TENS mini (1800528)



Način rada: TENS

Medio TENS (1800527)



Način rada: TENS

Medio TENS eco (1800508)



Način rada: TENS

Medio INCONTI (1800507)



Način rada: Inkontinencija • TENS

Medio EMS TENS eco (1800505)



Način rada: TENS • EMS

Medio IF eco (1800506)



Način rada: TENS • EMS • IF • MIC

Medio TENS micro (1800530)



Način rada: TENS

Medio TENS wireless (1800531)



Način rada: TENS

Medio SONO mini

(1800533)

Način rada:
Ultrazvuk 1 MHz




Medio UniSONO mini

(1800534)

Način rada:
Ultrazvuk 1 MHz •
TENS



Različiti pribor i potrošni materijal:



1. Podesiva elektroda 50x50 mm (1501741)
Podesiva elektroda 50x90 mm (1501742)
2. Gumena elektroda 45x50 mm (1501208)
Gumena elektroda 50x70 mm (1501209)
Gumena elektroda 80x115 mm (1501223)
3. Vlažan jastuk za gumenu elektrodu 45x50 mm (1502221)
Vlažan jastuk za gumenu elektrodu 50x70 mm (1502222)
Vlažan jastuk za gumenu elektrodu 80x115 mm (1502223)
4. Gumeni remen za fiksiranje 1 m (1504111)
5. Dugme za gumeni remen (1504201)
6. Kablovi za elektrode

Naziv uređaja:	Medio TENS mini	Medio TENS	Medio TENS eco	Medio Inconti	Medio EMS TENS eco	Medio IF eco	Medio TENS micro	Medio TENS wireless	Medio SONO mini	Medio UniSONO mini
Šifra	1800528	1800527	1800508	1800507	1800505	1800506	1800530	1800531	1800533	1800534

Načini rada:

TENS	✓	✓	✓	✓	✓	✓	✓	✓		✓
EMS					✓	✓				
IF						✓				
MIC						✓				
Inkontinencija				✓						
Ultrazvuk 1 MHz									✓	✓

Osnovni paket:

Podesiva elektroda 40x40 mm (1501701)	2 kom	2 kom	2 kom		2 kom	2 kom				
Podesiva elektroda 50x50 mm (1501741)	*	*	*	2 kom	*	*				
Podesiva elektroda 85x45 mm - oval							2 kom			
Samolepilne elektrode 90x60 mm - oval								2 kom		
Podesiva elektroda 50x100 mm										1 kom
Ultrazvučni gel									1 kom	1 kom
Vaginalna sonda				1 kom						
Gumena elektroda 45x50 mm (1501208)	*	2 kom	*	*	2 kom	*				
Vlažan jastuk za gumenu elektrodu 45x50 mm (1502221)	*	2 kom	*	*	2 kom	*				
Gumeni remen za fiksiranje 1 m (1504111)	*	2 kom	*	*	2 kom	*				
Kablovi za elektrode	2 kom	2 kom	2 kom	2 kom	2 kom	2 kom				1 kom
Baterija	✓	✓	✓	✓	✓	✓				
Adapter			*		*	✓			✓	✓

Tehnički podaci:

Broj kanala:	2	2	2	2	2	2	1	2	1	1
Frekvencija:	2 - 150 Hz	2 - 150 Hz	1 - 150 Hz	3 - 100 Hz	1 - 150 Hz	1 - 150 Hz	2 - 125 Hz	5 - 125 Hz	1 MHz	1 MHz, 1-120 Hz
Izlazna struja:	0 - 80 mA	0 - 80 mA	0 - 105 mA	0 - 90 mA	0 - 105 mA	0 - 105 mA	0 - 90 mA	0 - 90 mA	/	/
Širina impulsa:	30 - 260 μ s	30 - 260 μ s	50 - 300 μ s	50 - 450 μ s	50 - 300 μ s	50 - 300 μ s	100 - 200 μ s	130 - 250 μ s	/	/
Oblik pulsa:	Asimetričan bifazni, kružni	Asimetričan bifazni, kružni	Asimetričan bifazni, kružni	Asimetričan bifazni, kružni	Asimetričan bifazni, kružni	Asimetričan bifazni, kružni	Asimetričan bifazni, kružni	Asimetričan bifazni, kružni	/	Asimetričan bifazni, kružni
Napajanje / baterija:	9 V	9 V	9 V	4 x AAA	9 V	9 V	CR2032	6 x AAA	230 V	230 V
Dimenzije [mm]:	95x65x23	95x65x23	114x65x23	109x68x27	114x65x23	114x65x23	166x43x11	103x56x25	202x49x70	209x53x89
Težina:	110 g	110 g	127 g	105 g	127 g	127 g	11 g	45 g	193 g	235 g
Klasifikacija IEC 60601-1:	Klasa I, Tip BF									
Klasifikacija MDD 93/42/EEC:	Klasa IIb									

* Jedna zvjezdica označava dodatnu opremu

DISTRIBUTER:

PROIZVOĐAČ:

Iskra Medical d.o.o., Stegne 23, 1000 Ljubljana

Office: Otoče 5a, 4244 Podnart, Slovenia

T: +386 4 53 25 760, +386 4 53 25 762

F: +386 4 53 25 718, M: +386 41 622 918

E: info@iskramedical.eu, W: www.iskramedical.eu

